

O3elite

DUAL OZONE GENERATOR

Instructions: Find the ozone output you want in the chart below. Match the Oxygen LPM on the left with the O3Elite Dial Setting on top. For instance, for 5 gamma, set your oxygen to 3/4 LPM and your O3Elite Dial to 3. For 48 gamma, set your oxygen to 1/4 LPM, your O3Elite Dial to 5, and press Boost.

O2 LPM	DIAL SETTING							
	LOW	2	3	4	5	6	7	8
1	2	3	4	5	6	7	8	9
3/4	3	4	5	6	8	9	10	11
1/2	4	5	6	8	12	15	16	17
1/4	7	8	10	14	20	25	28	29
1/8	16	17	21	26	36	41	44	46
1/16	29	31	35	41	54	60	64	70

O2 LPM	DIAL SETTING + BOOST							
	LOW	2	3	4	5	6	7	8
1	4	6	9	12	15	16	17	18
3/4	7	8	11	14	18	20	21	22
1/2	11	12	16	22	28	31	32	33
1/4	22	24	29	38	48	50	52	53
1/8	35	38	45	54	61	66	70	74
1/16	65	67	70	77	84	88	93	100