

O3elite

SINGLE OZONE GENERATOR

Instructions: Find the ozone output you want in the chart. Match the Oxygen LPM on the left with the O3Elite Dial Setting on top. For instance, for 14 gamma, set your oxygen to 1/4 LPM and your O3Elite Dial to 4.

| O2 LPM | DIAL SETTING | | | | | | | |
|-----------|--------------|----|----|----|----|----|----|----|
| | LOW | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
| 1 | 2 | 3 | 4 | 5 | 6 | 7 | 8 | 9 |
| 3/4 | 3 | 4 | 5 | 6 | 8 | 9 | 10 | 11 |
| 1/2 | 4 | 5 | 6 | 8 | 12 | 15 | 16 | 17 |
| 1/4 | 7 | 8 | 10 | 14 | 20 | 25 | 28 | 29 |
| 1/8 | 16 | 17 | 21 | 26 | 36 | 41 | 44 | 46 |
| 1/16 | 29 | 31 | 35 | 41 | 54 | 60 | 64 | 70 |